

Anti-KIAA0101 antibody (Internal) (STJ72363)

STJ72363

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Goat polyclonal antibody anti-KIAA0101 (Internal) is suitable for use in ELISA and Western Blot research applications. |
| Applications | Pep-ELISA, WB-Trf |
| Host/Source | Goat |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Dilution Range | ELISA-antibody detection limit dilution 1:128000. |
| Formulation | 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA |
| Isotype | IgG |
| Storage Instruction | Store at -20 on receipt and minimise freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 9768 |
| Gene Symbol | PCLAF |
| Uniprot ID | PAF15_HUMAN |
| Immunogen | |
| Immunogen Region | Internal |
| Specificity | This antibody is expected to recognize isoform 1 (NP_055551.1) only. |
| Immunogen Sequence | SRKAENKYAGGN |