

Anti-Zfp652-Mouse antibody (C-Term) (STJ72689)

STJ72689

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Zfp652-Mouse (C-Term) is suitable for use in ELISA, Western Blot and Immunofluorescence/Immunocytochemistry research applications.
Applications	Pep-ELISA, WB, IF/ICC
Host/Source	Goat
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:16000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen Region	C-Term
Specificity	Reported variants represent identical protein: NP_055712.1, NP_001138837.1
Immunogen Sequence	RHLAEKNSAAPAQHH