

Anti-PCBP3 antibody (14-28) (STJ72726)

STJ72726

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-PCBP3 (14-28) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:32000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54039
Gene Symbol	PCBP3
Uniprot ID	PCBP3_HUMAN
Immunogen Region	14-28
Specificity	This antibody is expected to recognize both reported isoforms (NP_065389.2; NP_001123613.1). No cross-reactivity with other PCBP proteins are expected.
Immunogen Sequence	PHSTLSTLSHHPQPQ