

Anti-ZNF697 antibody (145-157) (STJ72923)

STJ72923

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-ZNF697 (145-157) is suitable for use in ELISA, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Pep-ELISA, WB-Trf, IHC, IP
Host/Source	Goat
Reactivity	Human, Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:128000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage	Store at -20 on receipt and minimise freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	90874
Gene Symbol	ZNF697
Uniprot ID	ZNF697_HUMAN
Immunogen	
Immunogen Region	145-157
Specificity	
Immunogen Sequence	RHLSLGSRRHGDK