

Anti-Centromere protein E antibody (2593-2604) (STJ73139)

STJ73139

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Centromere protein E (2593-2604) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:32000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1062
Gene Symbol	CENPE
Uniprot ID	CENPE_HUMAN
Immunogen	
Immunogen Region	2593-2604
Specificity	This antibody is expected to recognize both reported isoforms (NP_001804.2; NP_001273663.1).
Immunogen Sequence	KREAHKQVTCEN