

Anti-HSULF-1/sulfatase 1 antibody (Internal) (STJ73257)

STJ73257

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-HSULF-1/sulfatase 1 (Internal) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA, WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	ELISA-antibody detection limit dilution 1:128000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20 on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23213
Gene Symbol	SULF1
Uniprot ID	SULF1_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	Reported variants NP_055985.2; NP_001121676.1; NP_001121677.1; NP_001121678.1 represent identical protein
Immunogen Sequence	RGFHDKDKECS