

## Anti-Sgta-Mouse antibody (84-97) (STJ73523)

STJ73523

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Goat polyclonal antibody anti-Sgta-Mouse (84-97) is suitable for use in ELISA, Western Blot and Immunohistochemistry research applications.
<b>Applications</b>	Pep-ELISA, WB, IHC
<b>Host/Source</b>	Goat
<b>Reactivity</b>	Human, Mouse, Rat, Dog, Pig, Cow

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Dilution Range</b>	WB-1-2µg/ml IHC-5µg/ml ELISA-antibody detection limit dilution 1:2000.
<b>Formulation</b>	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20 on receipt and minimise freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	6449
<b>Gene Symbol</b>	SGTA
<b>Uniprot ID</b>	SGTA_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	84-97
<b>Specificity</b>	
<b>Immunogen Sequence</b>	PSEEDSAEAERLKT



STJ73523 (1µg/ml) staining of PD19 (A) and (2µg/ml) K562 (B) cell lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.