

## Anti-Phospho-Tubulin Alpha 1/3/4-Tyr272 antibody (210-290) (STJ90662)

STJ90662

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-Tubulin alpha-1A chain and Tubulin alpha-1B chain and Tubulin alpha-1C chain and OBSOLETE and Tubulin alpha-3E chain and Tubulin alpha-4A chain-Tyr272 (210-290) is suitable for use in Western Blot and ELISA
<b>Applications</b>	WB, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution</b>	WB 1:500-1:2000
<b>Range</b>	ELISA 1:40000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	7846 7277 112714 TUBA1A TUBA4A TBA1A_HUMAN TBA4A_HUMAN TBA3E_HUMAN
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human TUBA1/3/4 around the phosphorylation site of Tyr272 at amino acid range 238-287
<b>Immunogen Region</b>	210-290
<b>Specificity</b>	Phospho-Tubulin Alpha 1/3/4-Tyr272 polyclonal antibody (Tubulin alpha-1A chain and Tubulin alpha-1B chain and Tubulin alpha-1C chain and OBSOLETE and Tubulin alpha-3E chain and Tubulin alpha-4A chain) binds to endogenous Tubulin alpha-1A chain and Tu
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081