

## Anti-Phospho-VAV1-Tyr160 antibody (100-180) (STJ91129)

STJ91129

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Phospho-Proto-Oncogene Vav-Tyr160 (100-180) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications. |
| <b>Applications</b>      | IHC-P, IF-P, ELISA   |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human, Mouse, Rat  |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Polyclonal   |
| <b>Clone ID</b>            |  |
| <b>Concentration</b>       | 1 mg/mL  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.          |
| <b>Dilution Range</b>      | IHC 1:100-1:300<br>ELISA 1:5000  |
| <b>Formulation</b>         | PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.  |
| <b>Isotype</b>             | IgG  |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | <a href="#">7409</a>   |
| <b>Gene Symbol</b>        | <a href="#">VAV1</a>   |
| <b>Uniprot ID</b>         | <a href="#">VAV_HUMAN</a>  |
| <b>Immunogen</b>          | The antiserum was produced against synthesized peptide derived from human VAV1 around the phosphorylation site of Tyr160 at amino acid range 131-180                     |
| <b>Immunogen Region</b>   | 100-180  |
| <b>Specificity</b>        | Phospho-VAV1-Tyr160 polyclonal antibody (Proto-Oncogene Vav) binds to endogenous Proto-Oncogene Vav at the amino acid region 100-180 only when phosphorylated at Tyr160. |
| <b>Immunogen Sequence</b> |  |