

## Anti-Phospho-MARK1/2/3/4-Thr215 antibody (150-230) (STJ91168)

STJ91168

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Phospho-Serine/threonine-protein kinase MARK1 and Serine/threonine-protein kinase MARK2 and MAP/microtubule affinity-regulating kinase 3 and MAP/microtubule affinity-regulating kinase 4-Thr215 (150-230) is suitable for
<b>Applications</b>	IHC-P, IF-P, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution</b>	IHC 1:100-1:300
<b>Range</b>	ELISA 1:5000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">4139</a> <a href="#">57787</a> <a href="#">2011</a> <a href="#">MARK1</a> <a href="#">MARK4</a> <a href="#">MARK1_HUMAN</a> <a href="#">MARK4_HUMAN</a> <a href="#">MARK2_HUMAN</a>
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human MARK1/2/3/4 around the phosphorylation site of Thr215 at amino acid range 181-230
<b>Immunogen Region</b>	150-230
<b>Specificity</b>	Phospho-MARK1/2/3/4-Thr215 polyclonal antibody (Serine/threonine-protein kinase MARK1 and Serine/threonine-protein kinase MARK2 and MAP/microtubule affinity-regulating kinase 3 and MAP/microtubule affinity-regulating kinase 4) binds to endogenous Ser
<b>Immunogen Sequence</b>	