

## Anti-Actin-Alpha/Gamma antibody (40-120) (STJ91469)

STJ91469

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Actin, alpha cardiac muscle 1 and Actin, cytoplasmic 2 and Actin, gamma-enteric smooth muscle and Actin, alpha skeletal muscle (40-120) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence, Imm WB, IHC-P, IF, ICC, ELISA
<b>Applications</b>	WB, IHC-P, IF, ICC, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution</b>	WB 1:500-1:2000
<b>Range</b>	IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:1000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

#### Gene ID

72  
71  
58  
ACTG2  
ACTG1  
ACTH\_HUMAN  
ACTG\_HUMAN  
ACTS\_HUMAN  
<


**Immunogen** The antiserum was produced against synthesized peptide derived from human Actin at amino acid range 21-70

**Immunogen** 40-120


#### Region

**Specificity** Actin-Alpha/Gamma polyclonal antibody (Actin, alpha cardiac muscle 1 and Actin, cytoplasmic 2 and Actin, gamma-enteric smooth muscle and Actin, alpha skeletal muscle) binds to endogenous Actin, alpha cardiac muscle 1 and Actin, cytoplasmic 2 and Acti


#### Immunogen Sequence



Immunohistochemistry analysis of paraffin-embedded human skeletal muscle tissue, using Actin-pan Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from mouse brain, using Actin-pan Antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of HUVEC cells, using Actin-pan Antibody. The picture on the right is blocked with the synthesized peptide.