

Anti-ADAMTS17 antibody (150-230 Internal) (STJ91480)

STJ91480

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-A Disintegrin And Metalloproteinase With Thrombospondin Motifs 17 (150-230 Internal) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC-P, IF-P, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	IHC 1:100-1:300
Range	ELISA 1:20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	170691
Gene Symbol	ADAMTS17
Uniprot ID	ATS17_HUMAN
Immunogen	Synthesized peptide derived from ADAMTS-17. at amino acid range: 150-230
Immunogen Region	150-230 Internal
Specificity	ADAMTS17 polyclonal antibody (A Disintegrin And Metalloproteinase With Thrombospondin Motifs 17) binds to endogenous A Disintegrin And Metalloproteinase With Thrombospondin Motifs 17 at the amino acid region 150-230 Internal.
Immunogen Sequence	