

Anti-BRS3 antibody (140-220 Internal) (STJ91899)

STJ91899

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Bombesin Receptor Subtype-3 (140-220 Internal) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence, Immunocytochemistry and ELISA research applications.
Applications	WB, IHC-P, IF, ICC, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:10000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	680
Gene Symbol	BRS3
Uniprot ID	BRS3_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human BRS3 at amino acid range 161-210
Immunogen Region	140-220 Internal
Specificity	BRS3 polyclonal antibody (Bombesin Receptor Subtype-3) binds to endogenous Bombesin Receptor Subtype-3 at the amino acid region 140-220 Internal.
Immunogen Sequence	



Western blot analysis of lysates from HUVEC cells, using BRS3 Antibody. The lane on the right is blocked with the synthesized peptide.