

Anti-CNTNAP3 antibody (C-Term) (STJ92030)

STJ92030

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Contactin-Associated Protein-Like 3 (C-Term) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC-P, IF-P, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC 1:100-1:300 ELISA 1:10000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79937
Gene Symbol	CNTNAP3
Uniprot ID	CNTP3_HUMAN
Immunogen	Synthesized peptide derived from the C-terminal region of human CASPR3.
Immunogen Region	C-Term
Specificity	CNTNAP3 polyclonal antibody (Contactin-Associated Protein-Like 3) binds to endogenous Contactin-Associated Protein-Like 3 at the amino acid region C-Term.
Immunogen Sequence	