

## Anti-CD24 antibody (N-Term) (STJ92106)

## **GENERAL INFORMATION**

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-Signal Transducer Cd24 (N-Term) is suitable for use in Western Blot and ELISA research

applications.

Applications WB, ELISA Host/Source Rabbit

Reactivity Human, Mouse, Rat

## **PRODUCT PROPERTIES**

Clonality Polyclonal

Clone ID

Concentration 1 mg/mL Conjugation Unconjugated

**Purification** The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.

Dilution Range WB 1:500-1:2000

FLISA 1:10000

Formulation PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.

**Isotype** IgG

Storage Instruction Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

## **TARGET INFORMATION**

Gene ID 100133941

Gene Symbol CD24


Uniprot ID CD24\_HUMAN

Immunogen The antiserum was produced against synthesized peptide derived from human CD24 at amino acid range 14-63

Immunogen Region N-Term

Specificity CD24 polyclonal antibody (Signal Transducer Cd24) binds to endogenous Signal Transducer Cd24 at the amino acid region N-

Immunogen Sequence



-19 Western blot analysis of lysate from NIH/3T3 cells, using CD24 antibody.