

## Anti-CCP110 antibody (C-Term) (STJ92219)

STJ92219

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Centriolar Coiled-Coil Protein Of 110 Kda (C-Term) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
<b>Applications</b>	IHC-P, IF-P, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	IHC 1:100-1:300 ELISA 1:40000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">9738</a>
<b>Gene Symbol</b>	<a href="#">CCP110</a>
<b>Uniprot ID</b>	<a href="#">CP110_HUMAN</a>
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human CEP110 at amino acid range 850-900
<b>Immunogen Region</b>	C-Term
<b>Specificity</b>	CCP110 polyclonal antibody (Centriolar Coiled-Coil Protein Of 110 Kda) binds to endogenous Centriolar Coiled-Coil Protein Of 110 Kda at the amino acid region C-Term.
<b>Immunogen Sequence</b>	