

Anti-CEP41 antibody (Internal) (STJ92226)

STJ92226

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Centrosomal Protein Of 41 Kda (Internal) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC-P, IF-P, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC 1:100-1:300 ELISA 1:40000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	95681
Gene Symbol	CEP41
Uniprot ID	CEP41_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human CEP41 at amino acid range 150-200
Immunogen Region	Internal
Specificity	CEP41 polyclonal antibody (Centrosomal Protein Of 41 Kda) binds to endogenous Centrosomal Protein Of 41 Kda at the amino acid region Internal.
Immunogen Sequence	