

## Anti-CLDN3 antibody (140-220) (STJ92313)

STJ92313

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Claudin-3 (140-220) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications.
<b>Applications</b>	WB, IHC-P, IF-P, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES


<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:5000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


<b>Gene ID</b>	1365
<b>Gene Symbol</b>	CLDN3
<b>Uniprot ID</b>	CLD3_HUMAN
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human Claudin 3 at amino acid range 171-220
<b>Immunogen Region</b>	140-220
<b>Specificity</b>	CLDN3 polyclonal antibody (Claudin-3) binds to endogenous Claudin-3 at the amino acid region 140-220.
<b>Immunogen Sequence</b>	




Western blot analysis of lysates from COLO205 cells, treated with EGF 200ng/ml 30', using Claudin 3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of HeLa KB lysis using Claudin-3 antibody. Antibody was diluted at 1:2000



Western blot analysis of 293T 3T3 lysis using Claudin-3 antibody. Antibody was diluted at 1:2000



Western blot analysis of COLO cells using Claudin-3 Polyclonal Antibody diluted at 1: 2000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081