

## Anti-C1QTNF8 antibody (C-Term) (STJ92519)

STJ92519

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Complement C1q Tumor Necrosis Factor-Related Protein 8 (C-Term) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
<b>Applications</b>	IHC-P, IF-P, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution</b>	IHC 1:100-1:300
<b>Range</b>	ELISA 1:10000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">390664</a>
<b>Gene Symbol</b>	<a href="#">C1QTNF8</a>
<b>Uniprot ID</b>	<a href="#">C1QT8_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from the C-terminal region of human CTRP8.
<b>Immunogen Region</b>	C-Term
<b>Specificity</b>	C1QTNF8 polyclonal antibody (Complement C1q Tumor Necrosis Factor-Related Protein 8) binds to endogenous Complement C1q Tumor Necrosis Factor-Related Protein 8 at the amino acid region C-Term.
<b>Immunogen Sequence</b>	