

Anti-CYP4V2 antibody (Internal) (STJ92601)

STJ92601

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytochrome P450 4v2 (Internal) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	WB, IHC-P, IF-P, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	285440
Gene Symbol	CYP4V2
Uniprot ID	CP4V2_HUMAN
Immunogen	Synthesized peptide derived from the Internal region of human CYP4V2.
Immunogen Region	Internal
Specificity	CYP4V2 polyclonal antibody (Cytochrome P450 4v2) binds to endogenous Cytochrome P450 4v2 at the amino acid region Internal.
Immunogen Sequence	