

Anti-ITGB1 antibody (720-800) (STJ93733)

STJ93733

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Integrin Beta-1 (720-800) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-1:2000 ELISA 1:40000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3688
Gene Symbol	ITGB1
Uniprot ID	ITB1_HUMAN
Immunogen	Synthesized peptide derived from Integrin Beta 1. at amino acid range: 720-800
Immunogen Region	720-800
Specificity	ITGB1 polyclonal antibody (Integrin Beta-1) binds to endogenous Integrin Beta-1 at the amino acid region 720-800.
Immunogen Sequence	