

Anti-NFKB1 antibody (860-940) (STJ94462)

STJ94462

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nuclear Factor Nf-Kappa-B P105 Subunit (860-940) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC-P, IF-P, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC 1:100-1:300 ELISA 1:5000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4790
Gene Symbol	NFKB1
Uniprot ID	NFKB1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human NF-kappaB p105/p50 at amino acid range 891-940
Immunogen Region	860-940
Specificity	NFKB1 polyclonal antibody (Nuclear Factor Nf-Kappa-B P105 Subunit) binds to endogenous Nuclear Factor Nf-Kappa-B P105 Subunit at the amino acid region 860-940.
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded human ovary tissue, using NF-kappaB p105/p50 Antibody. The picture on the right is blocked with the synthesized peptide.