

## Anti-Olfactory receptor 3A2/3A3 antibody (40-120 Internal) (STJ94694)

STJ94694

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Olfactory receptor 3A3 and Olfactory receptor 3A2 (40-120 Internal) is suitable for use in Western Blot, Immunofluorescence, Immunocytochemistry and ELISA research applications.
<b>Applications</b>	WB, IF, ICC, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Rat, Mouse

### PRODUCT PROPERTIES


<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:5000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


<b>Gene ID</b>	<a href="#">4995</a> <a href="#">8392</a>
<b>Gene Symbol</b>	<a href="#">OR3A2</a> <a href="#">OR3A3</a>
<b>Uniprot ID</b>	<a href="#">OR3A2_HUMAN</a> <a href="#">OR3A3_HUMAN</a>
<b>Immunogen Region</b>	The antiserum was produced against synthesized peptide derived from human OR3A2/3 at amino acid range 65-114 40-120 Internal
<b>Specificity</b>	Olfactory receptor 3A2/3A3 polyclonal antibody (Olfactory receptor 3A3 and Olfactory receptor 3A2) binds to endogenous Olfactory receptor 3A3 and Olfactory receptor 3A2 at the amino acid region 40-120 Internal.
<b>Immunogen Sequence</b>	



Western blot analysis of lysates from MCF-7 cells, using OR3A2/3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using OR3A2/3 antibody.



Immunofluorescence analysis of MCF7 cells, using OR3A2/3 Antibody. The picture on the right is blocked with the synthesized peptide.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081