

## Anti-PER3 antibody (40-120 N-Term) (STJ95029)

STJ95029

### GENERAL INFORMATION


|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Period Circadian Protein Homolog 3 (40-120 N-Term) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications. |
| <b>Applications</b>      | WB, IHC-P, IF-P, ELISA  |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Human, Rat, Mouse   |

### PRODUCT PROPERTIES


|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Polyclonal   |
| <b>Clone ID</b>            |  |
| <b>Concentration</b>       | 1 mg/mL  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.          |
| <b>Dilution Range</b>      | WB 1:500-1:2000<br>IHC 1:100-1:300<br>ELISA 1:40000  |
| <b>Formulation</b>         | PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.  |
| <b>Isotype</b>             | IgG  |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION


|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | 8863   |
| <b>Gene Symbol</b>        | PER3   |
| <b>Uniprot ID</b>         | PER3_HUMAN   |
| <b>Immunogen</b>          | The antiserum was produced against synthesized peptide derived from human PER3 at amino acid range 21-70   |
| <b>Immunogen Region</b>   | 40-120 N-Term  |
| <b>Specificity</b>        | PER3 polyclonal antibody (Period Circadian Protein Homolog 3) binds to endogenous Period Circadian Protein Homolog 3 at the amino acid region 40-120 N-Term. |
| <b>Immunogen Sequence</b> |  |



Western blot analysis of lysates from Jurkat cells, treated with insulin 0.01U/ml 15', using PER3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using PER3 antibody.



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100 (4°C overnight). High-pressure and temperature Tris-EDTA, pH8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081