

Anti-PIAS2 antibody (10-90 N-Term) (STJ95084)

STJ95084

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E3 Sumo-Protein Ligase Pias2 (10-90 N-Term) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat, Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-1:2000 ELISA 1:20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	9063
Gene Symbol	PIAS2
Uniprot ID	PIAS2_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human PIAS2 at amino acid range 10-59
Immunogen Region	10-90 N-Term
Specificity	PIAS2 polyclonal antibody (E3 Sumo-Protein Ligase Pias2) binds to endogenous E3 Sumo-Protein Ligase Pias2 at the amino acid region 10-90 N-Term.
Immunogen Sequence	



Western blot analysis of COS7 cells using PIASx Polyclonal Antibody diluted at 1: 1000 cells nucleus extracted by Minute™ Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).



Western blot analysis of lysates from COS7 cells, using PIAS2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of various cells using PIASx Polyclonal Antibody diluted at 1: 1000 cells nucleus extracted by Minute™ Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081