

## Anti-PRRX1 antibody (30-110 Internal) (STJ95239)

STJ95239

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Paired Mesoderm Homeobox Protein 1 (30-110 Internal) is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES


<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:500-1:2000 ELISA 1:10000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">5396</a>
<b>Gene Symbol</b>	<a href="#">PRRX1</a>
<b>Uniprot ID</b>	<a href="#">PRRX1_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from PRRX1. at amino acid range: 30-110
<b>Immunogen Region</b>	30-110 Internal
<b>Specificity</b>	PRRX1 polyclonal antibody (Paired Mesoderm Homeobox Protein 1) binds to endogenous Paired Mesoderm Homeobox Protein 1 at the amino acid region 30-110 Internal.
<b>Immunogen Sequence</b>	



Western blot analysis of Jurkat cells using PRRX1 Polyclonal Antibody cells nucleus extracted by Minute™ Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).



Western blot analysis of various cells using PRRX1 Polyclonal Antibody cells nucleus extracted by Minute™ Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081