

## Anti-RAB24 antibody (90-170 Internal) (STJ95288)

STJ95288

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Ras-Related Protein Rab-24 (90-170 Internal) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
<b>Applications</b>	IHC-P, IF-P, ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
<b>Dilution Range</b>	WB 1:500-2000 IHC 1:100-1:300 ELISA 1:10000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	53917
<b>Gene Symbol</b>	RAB24
<b>Uniprot ID</b>	RAB24_HUMAN
<b>Immunogen</b>	Synthesized peptide derived from Rab 24. at amino acid range: 90-170
<b>Immunogen Region</b>	90-170 Internal
<b>Specificity</b>	RAB24 polyclonal antibody (Ras-Related Protein Rab-24) binds to endogenous Ras-Related Protein Rab-24 at the amino acid region 90-170 Internal.
<b>Immunogen Sequence</b>	