

Anti-UBE2G1 antibody (Internal) (STJ96164)

STJ96164

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ubiquitin-Conjugating Enzyme E2 G1 (Internal) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC-P, IF-P, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC 1:100-1:300 ELISA 1:40000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7326
Gene Symbol	UBE2G1
Uniprot ID	UB2G1_HUMAN
Immunogen	Synthesized peptide derived from the Internal region of human UBE2G1.
Immunogen Region	Internal
Specificity	UBE2G1 polyclonal antibody (Ubiquitin-Conjugating Enzyme E2 G1) binds to endogenous Ubiquitin-Conjugating Enzyme E2 G1 at the amino acid region Internal.
Immunogen Sequence	