

Anti-IFNA1/IFNA13 antibody (C-Term) (STJ96535)

STJ96535

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Interferon Alpha-1/13 (C-Term) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-1:2000 ELISA 1:20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3439/3447
Gene Symbol	IFNA1,IFNA13
Uniprot ID	IFNA1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human IFNA1/13 at amino acid range 140-189
Immunogen Region	C-Term
Specificity	IFNA1/IFNA13 polyclonal antibody (Interferon Alpha-1/13) binds to endogenous Interferon Alpha-1/13 at the amino acid region C-Term.
Immunogen Sequence	



Western blot analysis of BT474, A549 cells using IFN-Alpha 1 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081