

Anti-Ub antibody (STJ96797)

STJ96797

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ubiquitin-60S ribosomal protein L40 and Ubiquitin-40S ribosomal protein S27a and Polyubiquitin-B and Polyubiquitin-C is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-1:2000
Range	ELISA 1:20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID [6233](#)
[7316](#)
[7314](#)
[RPS27A](#)
[UBC](#)
[RS27A_HUMAN](#)
[UBC_HUMAN](#)
[UBB_HUMAN](#)

Immunogen Synthesized peptide derived from human Ub around the non-acetylation site of K63.

Immunogen Region

Specificity Ub polyclonal antibody (Ubiquitin-60S ribosomal protein L40 and Ubiquitin-40S ribosomal protein S27a and Polyubiquitin-B and Polyubiquitin-C) binds to endogenous Ubiquitin-60S ribosomal protein L40 and Ubiquitin-40S ribosomal protein S27a and Polyubi

Immunogen Sequence



Western blot analysis of HepG2 cells using Ub Polyclonal Antibody. Secondary antibody was diluted at 1:20000