

Anti-NTS antibody (121-170 C-Term) (STJ96900)

STJ96900

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Neurotensin/Neuromedin N Cleaved Into-Large Neuromedin N (121-170 C-Term) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-1:2000
Range	ELISA 1:20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4922
Gene Symbol	NTS
Uniprot ID	NEUT_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human NTS at amino acid range 121-170
Immunogen Region	121-170 C-Term
Specificity	NTS polyclonal antibody (Neurotensin/Neuromedin N Cleaved Into-Large Neuromedin N) binds to endogenous Neurotensin/Neuromedin N Cleaved Into-Large Neuromedin N at the amino acid region 121-170 C-Term.
Immunogen Sequence	



Western blot analysis of HeLa, HepG2 cells using Neurotensin Polyclonal Antibody. Secondary antibody was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081