

Anti-PTGS2 antibody [3D4] (STJ97061)

STJ97061

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Prostaglandin G/H Synthase 2 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Mouse
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3D4
Concentration	
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200
Formulation	PBS, pH 7.4, 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5743
Gene Symbol	PTGS2
Uniprot ID	PGH2_HUMAN
Immunogen	Synthetic peptide of COX2
Immunogen Region	
Specificity	PTGS2 monoclonal antibody (Prostaglandin G/H Synthase 2) binds to endogenous Prostaglandin G/H Synthase 2.
Immunogen Sequence	