

Anti-ECFP-Tag antibody (STJ97093)

STJ97093

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse polyclonal antibody anti-ECFP-Tag-Tag is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Tag

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:2000-5000
Formulation	PBS, pH 7.4, 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant Protein of ECFP-Tag
Immunogen Region	
Specificity	The antibody detects ECFP & CFP tag protein.
Immunogen Sequence	



Western blot analysis of ECFP Recombinant protein, diluted at 1) 1:5000, 2) 1:20000. Secondary antibody was diluted at 1:20000