

Anti-Acetyl-Ub-Lys33 antibody (STJ97363)

STJ97363

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-Ubiquitin-60S ribosomal protein L40 and Ubiquitin-40S ribosomal protein S27a and Polyubiquitin-B and Polyubiquitin-C-Lys33 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution	WB 1:500-1:2000
Range	ELISA 1:20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID [7314](#)
[6233](#)
[7316](#)
[UBB](#)
[RPS27A](#)
[UBB_HUMAN](#)
[RS27A_HUMAN](#)
[UBC_HUMAN](#)

Immunogen Synthesized acetyl-peptide derived from the human Ub around the acetylation site of K33.


Immunogen Region

Specificity Acetyl-Ub (K33) Polyclonal AntibodySynthesized peptide derived from the human Ub around the acetylation site of K33.

Immunogen Sequence



Western blot analysis of HeLa cells using Acetyl-Ub (K33) Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Western blot analysis of HeLa cells using Acetyl-Ub (K33) Polyclonal Antibody. Secondary antibody was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081