

Anti-PGLYRP1 antibody (131-180 Internal) (STJ97687)

STJ97687

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Peptidoglycan Recognition Protein 1 (131-180 Internal) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8993
Gene Symbol	PGLYRP1
Uniprot ID	PGRP1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human PGLYRP1 at amino acid range 131-180
Immunogen Region	131-180 Internal
Specificity	PGLYRP1 polyclonal antibody (Peptidoglycan Recognition Protein 1) binds to endogenous Peptidoglycan Recognition Protein 1 at the amino acid region 131-180 Internal.
Immunogen Sequence	



Western blot analysis of HepG2 823-AV 293T AD293 HeLa lysis using PGLYRP1 antibody. Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081