

Anti-FABP4 antibody (61-121) [5H11A4E11] (STJ97815)

STJ97815

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Fatty Acid-Binding Protein-Adipocyte (61-121) is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5H11A4E11
Concentration	
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Ascitic fluid, 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2167
Gene Symbol	FABP4
Uniprot ID	FABP4_HUMAN
Immunogen	Purified recombinant fragment of A-FABP (aa61-121) expressed in E.coli.
Immunogen Region	61-121
Specificity	FABP4 monoclonal antibody (Fatty Acid-Binding Protein-Adipocyte) binds to endogenous Fatty Acid-Binding Protein-Adipocyte at the amino acid region 61-121.
Immunogen Sequence	



Western blot analysis using A-FABP monoclonal antibody against truncated Trx-FABP4 recombinant protein (1).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081