

Anti-BID antibody [3C5] (STJ97864)

STJ97864

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Bh3-Interacting Domain Death Agonist is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence, Immunocytochemistry, Flow Cytometry and ELISA research applications.
Applications	WB, IHC-P, IF, ICC, FC, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	3C5
Concentration	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:500-1:2000 IHC 1:200-1:1000 IF 1:200-1:1000 FC 1:200-1:400 ELISA 1:10000
Formulation	Ascitic fluid, 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	637
Gene Symbol	BID
Uniprot ID	BID_HUMAN
Immunogen	Purified recombinant fragment of human BID expressed in E.coli.
Region	
Specificity	BID monoclonal antibody (Bh3-Interacting Domain Death Agonist) binds to endogenous Bh3-Interacting Domain Death Agonist.
Immunogen Sequence	



Flow cytometric analysis of HeLa cells using BID monoclonal antibody (green) and negative control (purple)



Immunohistochemistry analysis of paraffin-embedded prostate tissues (left) and tonsil tissues (right) with DAB staining using BID monoclonal antibody.



Immunofluorescence analysis of HeLa cells using BID monoclonal antibody (green), Blue: DRAQ5 fluorescent DNA dye, Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081