

## Anti-GATA5 antibody [1B11] (STJ98102)

STJ98102

### GENERAL INFORMATION



<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Transcription Factor Gata-5 is suitable for use in Western Blot and ELISA research applications.
<b>Applications</b>	WB, ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	1B11
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	WB 1:500-1:2000 ELISA 1:10000
<b>Formulation</b>	Ascitic fluid, 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	140628
<b>Gene Symbol</b>	GATA5
<b>Uniprot ID</b>	GATA5_HUMAN
<b>Immunogen</b>	Purified recombinant fragment of human GATA-5 expressed in E.coli.
<b>Immunogen Region</b>	
<b>Specificity</b>	GATA5 monoclonal antibody (Transcription Factor Gata-5) binds to endogenous Transcription Factor Gata-5.
<b>Immunogen Sequence</b>	



Western blot analysis using GATA-5 monoclonal antibody against HEK293 (1) and GATA5-hlgGfC transfected HEK293 (2) cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081