

Anti-PAX6 antibody [1C8] (STJ98310)

STJ98310

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Paired Box Protein Pax-6 is suitable for use in Western Blot, Flow Cytometry and ELISA research applications.
Applications	WB, FC, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1C8
Concentration	
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:500-1:2000 FC 1:200-1:400 ELISA 1:10000
Formulation	Ascitic fluid, 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5080
Gene Symbol	PAX6
Uniprot ID	PAX6_HUMAN
Immunogen	Purified recombinant fragment of human Pax-6 expressed in E.coli.
Immunogen Region	
Specificity	PAX6 monoclonal antibody (Paired Box Protein Pax-6) binds to endogenous Paired Box Protein Pax-6.
Immunogen Sequence	