

Anti-TYK2 antibody [8G8E9/8G9H5/8G8B3] (STJ98439)

STJ98439

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Non-Receptor Tyrosine-Protein Kinase Tyk2 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	8G8E9/8G9H5/8G8B3
Concentration	
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Ascitic fluid, 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7297
Gene Symbol	TYK2
Uniprot ID	TYK2_HUMAN
Immunogen	Purified recombinant fragment of Tyk 2 expressed in E.coli.
Immunogen Region	
Specificity	TYK2 monoclonal antibody (Non-Receptor Tyrosine-Protein Kinase Tyk2) binds to endogenous Non-Receptor Tyrosine-Protein Kinase Tyk2.
Immunogen Sequence	



Western blot analysis using Tyk 2 monoclonal antibody against truncated TYK2 recombinant protein (1) and Jurkat cell lysate (2).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081