

## Anti-ARRB1 antibody (STJ98467)

STJ98467

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Beta-Arrestin-1 is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	WB 1:1000-1:2000
<b>Formulation</b>	Buffer, 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% Sodium Azide, 50% Glycerol.
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	408
<b>Gene Symbol</b>	ARRB1
<b>Uniprot ID</b>	ARRB1_HUMAN
<b>Immunogen</b>	Purified recombinant human Arrestin-Beta-1 protein fragments expressed in E.coli.
<b>Immunogen Region</b>	
<b>Specificity</b>	ARRB1 monoclonal antibody (Beta-Arrestin-1) binds to endogenous Beta-Arrestin-1.
<b>Immunogen Sequence</b>	



Western blot analysis using Arrestin-Beta-1 monoclonal antibody against PC12, Mouse brain, Raw264.7 cell lysate.