

## Anti-L1CAM antibody (N-Term) (STJ98519)

STJ98519

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Neural Cell Adhesion Molecule L1 (N-Term) is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat, Dog, Rabbit

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	WB 1:1000-1:2000
<b>Formulation</b>	Buffer, 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% Sodium Azide, 50% Glycerol.
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	3897
<b>Gene Symbol</b>	L1CAM
<b>Uniprot ID</b>	L1CAM_HUMAN
<b>Immunogen</b>	Purified recombinant human NCAM-L1 (N-terminus) protein fragments expressed in E.coli.
<b>Immunogen Region</b>	N-Term
<b>Specificity</b>	L1CAM monoclonal antibody (Neural Cell Adhesion Molecule L1) binds to endogenous Neural Cell Adhesion Molecule L1 at the amino acid region N-Term.
<b>Immunogen Sequence</b>	



Western blot analysis using NCAM-L1 monoclonal antibody against HepG2 cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081