

Anti-HSP90AA1 antibody (STJ98573)

STJ98573

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Heat Shock Protein Hsp 90-Alpha is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:500-2000 ELISA 1:10000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3320
Gene Symbol	HSP90AA1
Uniprot ID	HS90A_HUMAN
Immunogen	Synthetic peptide of HSP 90
Immunogen Region	
Specificity	HSP90AA1 polyclonal antibody (Heat Shock Protein Hsp 90-Alpha) binds to endogenous Heat Shock Protein Hsp 90-Alpha.
Immunogen Sequence	