

Anti-MAPK10 antibody (STJ98873)

STJ98873

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Mitogen-Activated Protein Kinase 10 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-300
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5602
Gene Symbol	MAPK10
Uniprot ID	MK10_HUMAN
Immunogen	Recombinant Protein of JNK3
Immunogen Region	
Specificity	MAPK10 polyclonal antibody (Mitogen-Activated Protein Kinase 10) binds to endogenous Mitogen-Activated Protein Kinase 10.
Immunogen Sequence	



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JNK3 Rabbit pAb diluted at 1:200

Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JNK3 Rabbit pAb diluted at 1:200

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081