

Anti-COL4A1 antibody (STJ98923)

STJ98923

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Collagen Alpha-1 (Iv Chain Cleaved Into-Arresten is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC-P, IF-P
Host/Source	Mouse
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC-P 1:50-300
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1282
Gene Symbol	COL4A1
Uniprot ID	CO4A1_HUMAN
Immunogen	Synthetic peptide of Collagen IV
Immunogen Region	
Specificity	COL4A1 monoclonal antibody (Collagen Alpha-1 (Iv Chain Cleaved Into-Arresten) binds to endogenous Collagen Alpha-1 (Iv Chain Cleaved Into-Arresten.
Immunogen Sequence	