

Anti-NAIP antibody (1191-1240) (STJ99015)

STJ99015

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Baculoviral lap Repeat-Containing Protein 1 (1191-1240) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC-P, IF-P, ELISA
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	IHC-P 1:50-200 ELISA 1:10000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4671
Gene Symbol	NAIP
Uniprot ID	BIRC1_HUMAN
Immunogen	Synthetic peptide from human protein at amino acid range: 1191-1240
Immunogen Region	1191-1240
Specificity	NAIP polyclonal antibody (Baculoviral lap Repeat-Containing Protein 1) binds to endogenous Baculoviral lap Repeat-Containing Protein 1 at the amino acid region 1191-1240.
Immunogen Sequence	



Immunohistochemical analysis of paraffin-embedded Human-placenta, antibody was diluted at 1:100

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081