

Mouse Anti-Monkey IgG antibody {Gold} (STJ99381)

STJ99381

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Gold-conjugated mouse monoclonal anti-Monkey IgG secondary antibody. For use in most research applications.
Applications	WB, IHC-P, IF, ICC, ELISA
Host/Source	Mouse
Reactivity	Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	4-5 mg/mL
Conjugation	Gold
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:5000-10000 IF 1:5000-10000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	