

Rabbit Anti-Pig IgG Heavy & Light Chain antibody {Biotin} (STJ99409)

STJ99409

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Biotin-conjugated rabbit polyclonal anti-Pig IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB, IHC-P, IF, ICC, ELISA
Host/Source	Rabbit
Reactivity	Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	4-5 mg/mL
Conjugation	Biotin
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:5000-10000 IF 1:5000-10000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	