

Goat Anti-Treeshrew IgG Heavy & Light Chain antibody {FITC} (STJ99488)

STJ99488

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	FITC-conjugated goat polyclonal anti-Treeshrew IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB, IHC-P, IF, ICC, ELISA
Host/Source	Goat
Reactivity	Treeshrew

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	4-5 mg/mL
Conjugation	FITC
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Dilution Range	WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:5000-10000 IF 1:5000-10000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	