

Goat Anti-Treeshrew IgG Heavy & Light Chain antibody {HRP} (STJ99538)

STJ99538

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Secondary antibodies |
| Short Description | HRP-conjugated goat polyclonal anti-Treeshrew IgG Heavy & Light Chain secondary antibody. For use in most research applications. |
| Applications | WB, IHC-P, IF, ICC, ELISA |
| Host/Source | Goat |
| Reactivity | Treeshrew |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 4-5 mg/mL |
| Conjugation | HRP |
| Purification | The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography. |
| Dilution Range | WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:5000-10000 IF 1:5000-10000 |
| Formulation | PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |